



IdentiPol QA2 application for biomedical material assessment

This abstract and poster paper show how the IdentiPol was used to measure the glass transition temperatures of the formulations in order to ensure they had the correct performance used in Body-temperature reverting, plasticised shape-memory polymers for the healing of fractures.

COMPRESSION SCREWS FOR THE FIXATION OF SMALL BONES AND BONE FRAGMENTS USING BIOCOMPATIBLE, BIORESORBABLE, BODY-TEMPERATURE REVERTING, PLASTICISED SHAPE MEMORY POLYMERS

B.M. Thomson¹, G. Thompson¹, D. Vgenopoulos¹, K. Nair¹, J. Duncan², K. Howell¹, M. Martyn¹ and P. Coates¹

¹School of Engineering and Informatics, Bradford University ²Lacerta Technology, Leicestershire

Orthopaedic compression screws promote fracture healing by drawing-together and stabilising adjacent bone fragments (e.g. Herbert screws in scaphoid repair).

Shape memory polymers are materials that switch from one macroscopic shape to another following the application of external 'triggering' stimuli (e.g. by shrinking in length when heated). A shape memory compression screw could therefore potentially be used to pull bone fragments together.

The shape memory phenomenon occurs because when polymers solidify from a melt (or dry from solutions), the molecules form tangled networks in which individual molecules adopt conformations with the lowest available energies. If a solid polymer is (i) heated to a temperature above the point where the polymer chains are able to begin to move relative to one another (called T_g) but below its melting point; (ii) stretched; and then (iii) rapidly cooled below T_g without releasing the tension, then the polymer's chains become 'locked' into energetically unfavourable conformations, (analogous to frozen, stretched rubber bands). As long as the polymer remains below T_g then its molecules remain locked in these energetically unfavourable conformations. However, if the polymer is heated above T_g, then molecules will move relative to one another and revert to their original energetically favourable conformations (e.g. re-heating a stretched shape memory polymer above T_g will cause it to shrink in length).

The biocompatible, bioresorbable polymer polylactic acid can show shape memory properties and could therefore potentially be used to make a shape memory compression screw. Unfortunately, the triggering temperature required to initiate shape reversion in unmodified PLA is too high to be conveniently incorporated into routine clinical practice.

The value for T_g for a particular polymer can be reduced by adding plasticisers, low molecular weight compounds that facilitate the movement of molecular chains within solid polymers. To identify plasticisers that can reduce the T_g of PLA to 37C (i.e. physiological temperatures), standardised mixtures of PLA, plasticiser and solvent were added to moulds and dried in a computer controlled oven to produce plasticiser-polymer films. The T_g of these materials was determined using an identiPol machine. The films' shape memory properties were assessed by (i) heating 10 mm long pieces of film to 85C in water; (ii) stretching them to 20mm long; (iii) rapidly cooling them to 5C without releasing the tension; and then (iv) attempting to trigger shape memory reversion by brief immersion in water at 10C – 60C or more prolonged immersion in saline at 37C.

Results showed that two reportedly biocompatible plasticisers, Triacetin and Tributyl *O*-acetylcitrate reduced PLA's T_g in a dose dependant manner and that shape memory reversion could occur at 37C in saline. Such materials could potentially be used to form compression screws that shrink in length post-implantation and pull bone fragments together.

Ongoing work aims to optimise the conditions required to produce the shape memory plasticised PLA. Future work will use the manufacturing and characterisation equipment available at Bradford (e.g. injection moulding, extrusion and machining) to produce compression screws from these shape memory materials.

COMPRESSION SCREWS FOR THE FIXATION OF SMALL BONE FRAGMENTS USING BIORESORBABLE, BODY-TEMPERATURE REVERTING, PLASTICISED SHAPE-MEMORY POLYMERS


B.M. Thomson¹, G. Thompson¹, D. Vgenopoulos¹, K. Nair¹, J. Duncan², K. Howell¹, M. Martyn¹
and P. Coates¹

Faculty of Engineering & Informatics, School of Engineering, University of Bradford¹, Lacerta Technology²

Introduction

OBJECTIVE: Make orthopaedic compression screws from stretched Shape Memory Polymers (SMP) and promote bone healing by using SMP's reversion force to gently pull bone fragments together

Underlying principle behind Shape Memory Polymers:



- The biocompatible, bioresorbable polymer polylactic acid (PLA) can show shape memory effects, but its T_g is too high for routine clinical use ($\sim 65^\circ\text{C}$).
 - Plasticisers are small molecules that reduce a polymer's T_g
- Find a biocompatible PLA-Plasticiser combination that reverts under physiological conditions and use this material to make an orthopaedic compression screw

Concept:

Stretched SMP screw on implantation

①

Post-implantation

②

Reversion above T_g




Scaphoid nonunion

→ SMP's reversion force gently pulls bone fragments together


Experimental

Identify a biocompatible Polylactate-Plasticiser SMP that reverts under physiological conditions

Screening method for measuring T_g of polymer-plasticiser-solvent cast films



Screening method for assessing reversion of polymer-plasticiser SMP films



For implants that exploit SMP phenomenon clinically, we want shape-memory polymers that are:

- Dimensionally stable at room temperature when dry
- Revert under physiological conditions (i.e. saline and 37°C)

Results

Design of compression screw



Typical existing design: Synthes 3.0mm Headless Compression screw, long thread and Solidworks model based upon several existing designs

Possible design criteria

THREAD PITCH




Thread pitch different at proximal and distal ends to produce compression

THREAD PROFILES



Symmetrical profile




Reverse buttress and buttress thread profiles




German Buttress Thread $33^\circ 3^\circ$


Use of biocompatible plasticiser to lower the T_g of polylactate so that SMP reversion occurs under physiological conditions




Plasticisers



N-Methyl-2-Pyrrolidone




Triacetin



Tributyl *O*-acetyl citrate (TAC)



Safety/ regulatory: Triacetin and Tributyl *O*-acetyl citrate are not hazardous according to Regulation (EC) No. 1272/2008 or classified as dangerous according to Directive 67/548/EEC.

Effect of plasticisers on T_g of polylactic acid (200mg)





Effect of Tributyl *O*-acetyl citrate (TAC) on reversion temperature of polylactate SMP

Reversion of plasticised polylactate SMP following 10 second immersion in water at 0-60°C



Reversion of plasticised polylactate SMP films (200mg PLA + NN μl TAC) at 37°C in saline



Conclusions

- Results showed that two reportedly biocompatible plasticisers, Triacetin and Tributyl *O*-acetyl citrate could reduce PLA's T_g in a dose dependant manner
- Shape memory polymers formed from Tributyl *O*-acetyl citrate plasticised PLA reverted under physiological conditions (37°C in saline).
- Such materials could potentially be used to form compression screws that shrink in length post-implantation and pull bone fragments together.